

GOL SITE

دلتا کیر
DELTA CARE

GPP-100

FDA for Analyzer

Automatic Specific Protein Analyzer



FEATURES

دلٲا كير
DELTA CARE

GOL SITE



All in one
Cartridge



One Step




Reliable
Result



Compact

- Classical Immunoassay Technology: Nephelometry
- All-in-one cartridge for single test
- Auto ID of assay QR code on reagent cartridge
- Small sample volume and wide test menu
- Quantitative measurement and maintenance free
- Independent calibration



OPERATION

ONE STEP ONLY!




TEST MENU

Assay	Linear Range	Sample Type	Category
HbA1c	2.8% – 15.5%	Blood	Diabetes
IgA	0.23 – 5.10 g/L	Blood/Plasma/Serum	Immune Disease
IgG	1.8 – 26.0 g/L	Blood/Plasma/Serum	
IgM	0.18 – 4.60 g/L	Plasma/Serum	
C3	0.19 – 4.70 g/L	Plasma/Serum	
C4	0.06 – 2.20 g/L	Plasma/Serum	
ASO	41 – 607 IU/mL	Serum	Rheumatoid Disease
RF	20 – 306 IU/mL	Serum	
Anti-CCP	15 – 200 IU/mL	Serum	
D-Dimer	0.28 – 15.70 mg/L	Plasma	Cardiovascular Disease
UsCRP	0.22 – 19.00 mg/L	Plasma/Serum	
LP(a)	60 – 1200 mg/L	Serum	
FR-CRP	0.5 – 380.0 mg/L	Blood/Plasma/Serum	Inflammation
CRP	1.6 – 152.0 mg/L	Plasma/Serum	
SAA	3 – 200 mg/L	Blood/Plasma/Serum	
Urine IgG	4 – 100 mg/L	Urine	Kidney Disease
Urine TRF	1 – 80 mg/L	Urine	
A1M	5 – 200 mg/L	Urine	
mALB	9.5 – 230.0 mg/L	Urine	
CysC	0.4 – 7.6 mg/L	Plasma/Serum	
BMG	0.4 – 13.0 mg/L	Plasma/Serum/Urine	


More assays are to be expected

CLINICAL PERFORMANCE

GPP-100 VS Biorad D10



GPP-100 VS Beckman IMMAGE



Excellent correlation result compared with well-established brand



SPECIFICATIONS

Measuring Time	3 to 8 minutes	Display	Color touch screen
Sample Volume	3 μ L to 300 μ L	Printer	Built-in thermal printer
Pack Size	20 T/Kit	Net Weight	7.1 kg
Memory Storage	Up to 100,000 test results	Size	404 mm*161 mm*397 mm

الوكيل الوحيد بجمهورية مصر العربية
شركة

For Medical Supplies
01010060298 - 01210629090
www.delta-care.com

دلتا كير
DELTA CARE